

Plexus Caudal Support Pads



Made in the U.S., Austin Edens' Plexus Pads feature a patent pending, common sense design.

Plexus Open Toe Caudal Support pads are designed to fit both fronts and hinds. The unique diamond-shaped holes provide greater locking strength for **FootPro™ DIM** and **Vettec Equi-Pak**. The extra length and width at the toe and heels accommodates any shaped foot and shoe modifications. The pads include a large frog plate, with heel check relief, for easy application on contracted heels. Frog plate stabilizers keep its tip anchored in place - even in the harshest conditions. The wedge pad features a 3-degree wedge for help in correcting negative palmar angles and tendon issues, while maintaining caudle support.



KERCKHAERT
FARRIER PRODUCTS

P.O. Box 1328
Shelbyville, KY 40066

www.farrierproducts.com

www.bellotafarrier.com



@farrierproductdistribution

@fpdinc

@farrierproducts



What's
New in Hoof
and Frog
Health




Common
Sense
Design

RELY ON THESE PRODUCTS FOR COMMON SENSE HOOFCARE SOLUTIONS

Watch videos of hoofcare
products in action at:
tinyurl.com/yc8mc3np



FootPro™
Cast,
Liberty
Cu Nails



FootPro™ Cast

- Resin impregnated fiberglass mesh fabric conforms to the hoof, provides support, wear protection and stabilizes the whole hoof.
- Dry application technique provides much greater working time and control by allowing the user the choice when to activate and cure the resin.
- Two sizes gives the user a choice depending on hoof size.



Allen's Blue
Powder,
FootPro™
DIM 20



FootPro™ DIM

- Non-toxic, non-allergenic formula that does not adversely affect structures of the hoof, nor is it affected by moisture.
- Easy to mix two equal parts by hand.
- Choose between cushion (DIM20) and support (DIM40), both available in 2lb. and 4lb. sets.



Plexus Pad,
FootPro™
DIM 20,
Kerckhaert
Aluminum
Triumph,
Liberty
Carrera
Nails



Plexus Pads

- Shape allows for easy fit for front or hind, reducing need for more inventory.
- Design also allows fit for more hoof sizes, also reducing need for more inventory.
- Frog support is tapered at end of pad to reduce the need to trim or grind pads, or grind heels of shoes.
- Open Toe design avoids trapping of foreign materials in front half of foot.
- High quality injection molded material should allow for resets.

